

As A Women, We Care For Our Ovaluation Period

Women's Guard U-TEST Saliva Ovulation Device For Smartphone is to observe the density of fern crystallization from your saliva(NaCl) to determine whether or not has ovulated. Multiple and serial uses should be used for this purpose.



Non-fertile Phase Pebbles

During this period of time, you will only see dot-shape crystals from your smartphone.



Fertile Phase Ferns

During this period of time, you will see fern-shape crystals covering the lens from your smartphone. This tells you that you are entering the ovulation period. Suggest to take photos every time you take testing and save them in your APP. They can be your reference for next ovulation period.



Fern-shape crystals

Transition Period

During this period of time, you will see some fern-shape crystals from your smartphone. This tells you that you are close to the peak of ovulation after 3 to 4 days. Our APP will remind you to use U-TEST for testing.

TARGET USERS

- Women who like to use a natural way to track their own menstrual period.
- Women who have irregular menstrual period.
- Middle age women before entering menopause who need a long period of time to track their hormone level.
- Women who prefer to use non-invasive or reusable device to do test.
- Women who don't want to use IUD(Intrauterine Device) or any other method to become pregnant.
- A couple who wants to increase the chance to have children.
- Women can do spa or facial treatment based on the correct period data.


FEATURES OF U-TEST

- Easy to Use
- Work with Any Smart Phones
- Free APP for Download
- Never Miss Your Test
- Non-Invasive Method
- A Trendy and Lightweight Device to Carry With
- Cost-Effective
- Scratch Resistance(Lens)
- Patents Granted
- U-Test Saliva Ovulation Device is non-invasive; it is the most natural method to test your fertile and non-fertile phase.
- U-Test Saliva Ovulation Device is granted with multiple patents worldwide, it also wins the award in Life Science Category of Red Dot Award in 2016.
- U-Test Saliva Ovulation Device is designed developed by the R&D team of Women's Guard Biomedical Technology Corp., Dr. Zheng, the outstanding Obstetrics and Gynecology Doctors of Kaohsiung Medical University is invited to participate during the development.



Use with U-Test APP


U-Test Saliva Ovulation Device can work with your smartphone to take pictures, store images, and set up your personal alarm as testing reminder. You can manage and oversee your menstrual period efficiently.




Press OK button when you see the crystal pattern after taking a picture.




Auto link to the period record board and write down the special code for the pattern.




A short video to guide how to use the device if this is your first time.



If you have any questions about the product, you can find the answers here.




The device has a reminder with all kind of information for your review.



Find out your fertile phase quickly by data.



What are the symptoms of ovulation? You can find the answers here.



Set up your personal alarm for testing reminding in Setting.

Patents Granted



Compatibility

Requires Android 4.4 or later

Requires IOS 8 or later
iphone 4s or later



Panpac Medical Corporation TEL : 886-2-86472242 www.panpac.com.tw
6F-2, No.202, Sec.3, Ta-Tong Road, Shi-Chih District, New Taipei City, Taiwan.